

MEDICAL CARE KIT FOR MEDICAL PERSONNEL

| Nr | Item | Description | Pieces |
|----|---------------------|---|--------|
| 1 | Biohazard Waste Bag | <ul style="list-style-type: none"> • Biohazard Waste Bag is made of yellow High Density Polyethylene. • Provides strength and durability with a capacity of 4-5 liters | 1 |
| 2 | Nitril gloves | <ul style="list-style-type: none"> • Is made of nitrile • size L • single use • CE mark. • expires after min 2 years from the manufacturing date. | 8 |
| 3 | Tourniquet | <ul style="list-style-type: none"> • Ensures the compression of the affected limb by twisting a rigid stick attached to the tourniquet. • Components: compression strap, rigid stick, safety straps • The compression strap a resistant material which generates pressure required for hemostasis when under tension. <ul style="list-style-type: none"> ○ It is not elastic. ○ Can be easily cleaned ○ it is water resistant • the rigid stick is made of a resistant material (metal or similar) which is attached to the tourniquet and when twisted generates tension in the compression strap. • the safety strap allows a good fixation of the rigid stick in order to maintain the pressure generated by the compression strap. It can be easily attached and easily detached • the tourniquet can be used one handed • the tourniquet is packed in water-resistant packaging, ready to use • single use | 2 |

| | | | |
|---|------------------------------|---|---|
| | | <ul style="list-style-type: none"> individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time,). CE mark. expires after min 2 years from the manufacturing date. | |
| 4 | Big Cinch Abdominal Bandaged | <ul style="list-style-type: none"> Elastic wrap made with natural rubber latex Dimensions: 30x40 cm single use individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). CE mark. expires after min 2 years from the manufacturing date. | 1 |
| 5 | Compressive dressing | <ul style="list-style-type: none"> Creates an effective compression by twisting a dressing around a clamp attached to the elastic dressing. Components: dressing, rigid clamp the dressing is made of an elastic material which allows a good compression of the area on which is used. the clamp is made of a resistant rigid material which allows to generate tension in the dressing and also maintains the dressing in position for the best compression single use expires after min 2 years from the manufacturing date. | 7 |
| 6 | KERLIX Gauze Bandage Rolls | <ul style="list-style-type: none"> is made of cotton Dimensions: 11 cm x 10 m single use individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). CE mark. expires after min 2 years from the manufacturing date. | 1 |
| 7 | Thoracic valve dressing | <ul style="list-style-type: none"> Used to drain the air from a tension pneumothorax. Adhesive dressing with a unidirectional valve. | 1 |

| | | | |
|----|------------------------|---|---|
| | | <ul style="list-style-type: none"> • single use • individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date. | |
| 8 | Oro-pharyngeal airway | <ul style="list-style-type: none"> • Specially designed tube that after insertion through the nose ensures the maintenance of an open airway in an unconscious patient. • length 90 mm and 100mm • made of latex free material which allows bending and smooth insertion in the nose • the outer end has a safety system which prevents the tube to be completely accidentally introduced or aspired • single use • CE mark. • expires after min 2 years from the manufacturing date. | 2 |
| 9 | Nasopharyngeal airway | <ul style="list-style-type: none"> • Specially designed tube that after insertion through the nose ensures the maintenance of an open airway in an unconscious patient. • 28 Fr diameter • length 165 - 175 mm • made of latex free material which allows bending and smooth insertion in the nose • the outer end has a safety system which prevents the tube to be completely accidentally introduced or aspired • single use • sterile • individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date. | 2 |
| 10 | TMM-CR, Cricothyrotomy | <ul style="list-style-type: none"> • content: - Tracheostomy Tube Cuffed, 6.0 mm – 1 pcs | 1 |

| | | | |
|----|----------------------|---|----|
| | | <ul style="list-style-type: none"> - Hemostat, Kelly Forceps Curved – 1 pcs - Scalpel, Sterile – 1 pcs - Syringe, Sterile, 10 mL – 1 pcs - Tracheal Hook, 4 - 1 pcs <ul style="list-style-type: none"> • expires after min 2 years from the manufacturing date. | |
| 11 | Ambu bag | <ul style="list-style-type: none"> • Ambu bags are respiration device used to provide positive pressure for manual ventilation needed when a patient lacks respiratory drive during resuscitation involving oxygen • Adult dimension • CE mark. • expires after min 2 years from the manufacturing date. | 1 |
| 12 | Oxygen mask | <ul style="list-style-type: none"> • Oxygen Mask is used for providing oxygen in the respiratory system through the mouth • Adults dimension | 1 |
| 13 | Decompression Needle | <ul style="list-style-type: none"> • sterile decompression needle designed for the removal of air from the pleural cavity • sterile • single use • individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date. | 2 |
| 14 | Burn Dressing | <ul style="list-style-type: none"> • dressing used to treat a variety of skin injuries including first, second, and third-degree burns whether thermal, chemical or electric. • 2 different dimensions: 10 cm x 40 cm and 10 cm x 10cm • sterile • single use • individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date. | 2 |
| 15 | Non- sterile | <ul style="list-style-type: none"> • material made of cotton 100%, with a good absorbent capacity | 10 |

| | | | |
|----|--------------------|--|---|
| | Dressings | <ul style="list-style-type: none"> • 6-8 sheets, • dimensions: 5 x 5 cm. • non-sterile • pairs are individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, etc.). • CE mark. • expires after min 2 years from the manufacturing date | |
| 16 | Adhesive tape | <ul style="list-style-type: none"> • Adhesive tape used to attach different gausses or dressings • for medical use • single use • individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date. | 2 |
| 17 | Triangular bandage | <ul style="list-style-type: none"> • Made by a soft but resistant material • used for basic splinting of the upper limbs • one angle of the equilateral triangle is 90⁰ and the two equal sides have 50-70 cm length | 2 |
| 18 | Medical Splint | <ul style="list-style-type: none"> • Splint is adaptable for use in almost any situation on almost any part of the body. • is built from a thin core of aluminum alloy, sandwiched between two layers of closed-cell foam • dimension: 90 cm • | 2 |
| 19 | IV Catheter | <ul style="list-style-type: none"> • IV Catheter small, flexible tube designed to deliver IV medications and fluids to the patient • Dimension 1.7 x 50 mm • Sterile • single use • individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date. | 4 |

| | | | |
|----|----------------------------|---|----|
| 20 | Eye plasters | <ul style="list-style-type: none"> • Small gauzes that allows protection of the wounded eyes • individually packed • single use • water resistant material on the external side and absorbent on the interior. • single use • individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date. | 4 |
| 21 | Antiseptic wipes | <ul style="list-style-type: none"> • Wipes with clorhexidine • packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date | 30 |
| 22 | Plasters – different sizes | <ul style="list-style-type: none"> • the plasters are padded. • single used • 5 different dimensions: the padded area will be at least 1 x 1 cm up to 10 x 10 cm • the adhesive area will ensure a good adhesion of the plaster to any part of the body and a good sealing of the affected area • single use • individually packed • single use • individually packed as regulated by EU-MDD (EU Medical Device Directive) (name of the product, description, expiration time, sterilization method, etc.). • CE mark. • expires after min 2 years from the manufacturing date. | 10 |
| 23 | Permanent Marker | | 1 |
| 24 | Scissors | <ul style="list-style-type: none"> • Bandage scissors are angled tip scissors, with a blunt tip on the bottom blade, used for sizing | 1 |

| | | | |
|----|-------|--|---|
| | | dressings and removing circumferential bandages | |
| 25 | Pouch | <ul style="list-style-type: none">• low-profile medical backpack.• compatible with PALS (Pouch Attachment Ladder System) webbing.• the pouch is water resistant• The strap and the pouch has a good mechanical resistance to tearing or friction• the color will be decided by the beneficiary and discussed with the provider | 1 |